



# Manual of ToTEM Software

**Prof. Fang Lin**

**South China Agricultural University**

**E-mail: [flin\\_rew@163.com](mailto:flin_rew@163.com)**

# CONTENT

---

- I. QUICK OPERATION.
- II. ENVIRONMENT CONFIGURATION.
- III. DETAILS






# I. QUICK OPERATION


0. Run `mex_cuda.m` and then `ToTEM_submit_v1.m` (source code).

1、 Read information about a supercell.



Anyone of three

2、 Press radio button to get a special direction.




Anyone of two

3a、 Slicing crystal by pressing 'Slicing' radio button to control a sub-dialog.




3b. Set the total slice number and pick some slices.





#### 4、 Selection for atomic scattering factor and vibration. For example:



The screenshot shows a software window with the following elements:



- Projected Potential** section:
  - Potential param.**: A dropdown menu with a blue arrow pointing to it from the text "select atomic scattering factor". The menu is open, showing four options: "Peng 0~6 A-1", "Peng 0~6 A-1", "Peng extended" (highlighted in blue), and "Lobato".
  - Vibration**: A radio button.
  - Temperature Factor**: A radio button.
  - Height Influence in Slice**: A checkbox.
- Imaging Region**: A section at the bottom.

Other visible elements include a value of "30" next to the "Vibration" radio button, a value of "2" next to the "Height Influence in Slice" checkbox, and a blue bracket on the right side of the "Vibration" and "Temperature Factor" options.

select atomic scattering factor

Choose whether to consider vibration or temperature factor, whether to assign potential into multiple slices.

5a、 Select STEM/HRTEM mode and set region and parameters. A STEM example:



Imaging Condition

Voltage: 300 kV Convergence: 10 mrad

Focus Spread: 3.0 nm Number: 1

Spatial Convergence: 0.3 mrad

Beam tilt: 0 mrad 0 degree

Magnitude(nm) & Angle(degree)

Defocus C1: -58.5 nm

2-fold astigmatism A1: 0 nm 0

3-fold astigmatism A2: 0 nm 0

Axial coma B2: 0 nm 0

(C3) Cs: 1200 um

4-fold astigmatism A3: 0 um 0

Axial star aberration S3: 0 um 0

5-fold astigmatism A4: 0 um 0

3-lobe aberration D4: 0 um 0

coma B4: 0 um 0

C5: 0 mm

6-fold astigmatism A5: 0 mm 0

R5: 0 mm 0

S5: 0 mm 0

STEM Detector mrad

40,200

5. Simulation

Path Save as: test



Other Parameter Settings and simulation.



5b、 Select STEM/HRTEM mode and set region and parameters. A HRTEM example.

Imaging Region

Edge Size  Pixels **4. Show**

Imaging Resolution  A/Pixel

Imaging Region (Top-Left)

Angstrom


Width & Height of Imaging Region

Pixels

3. Mode

☒ HRTEM ☐ CBED ☐ CPU

☐ STEM ☐ IDPC ☒ GPU



Imaging Condition

Voltage:  kV Aperture  1/nm

Focus Spread:  nm Number:

Spatial Convergence:  mrad

Beam tilt:  mrad  degree

Magnitude(nm) & Angle(degree)

Defocus C1:  nm

2-fold astigmatism A1:  nm

3-fold astigmatism A2:  nm

Axial coma B2:  nm

(C3) Cs:  um

4-fold astigmatism A3:  um

Axial star aberration S3:  um

5-fold astigmatism A4:  um

3-lobe aberration D4:  um

coma B4:  um

C5:  mm

6-fold astigmatism A5:  mm

R5:  mm

S5:  mm

**5. Simulation**

Path Save as:



## II. ENVIRONMENT CONFIGURATION

### 1. Installer


You can run this software on Windows or Linux with or without a CUDA environment. Here's how to run it on Windows and Linux (Ubuntu) with a CUDA environment.

Run `mex_cuda.m` first!

#### 1.1 WINDOWS

##### 1.1.1 CUDA


Open <https://developer.nvidia.com/cuda-downloads>  
Installation is successful if display as Fig. 1-1,



```
Windows PowerShell
版权所有 (C) Microsoft Corporation。保留所有权利。
尝试新的跨平台 PowerShell https://aka.ms/pscore6

PS C:\Users\YUANPENGJUN> nvcc -V
nvcc: NVIDIA (R) Cuda compiler driver
Copyright (c) 2005-2020 NVIDIA Corporation
Built on Mon_Oct_12_20:54:10_Pacific_Daylight_Time_2020
Cuda compilation tools, release 11.1, V11.1.105
Build cuda_11.1.relgpu_drvr455TC455_06.29190527_0
PS C:\Users\YUANPENGJUN> |
```

Type 'nvcc-v' at the terminal.



### 1.1.2 MCR

Open <https://www.mathworks.com/products/compiler/matlab-runtime.html> to install the 'MATLAB Runtime Version 2020b (9.9) Version'. Please choose the version on 'Windows 64-bit'.

### 1.1.3 RUN

Double-click 'ToTEM.exe' to run it.

## 1.2 Ubuntu

### 1.2.1 CUDA

Open <https://developer.nvidia.com/cuda-downloads> to install 'cuda 11 version', Then type 'nvcc-v' at the terminal. If the terminal has CUDA 11 version information, installation is successful.

### 1.2.2 MCR

Open <https://www.mathworks.com/products/compiler/matlab-runtime.html> to install the 'MATLAB Runtime Version R2019b (9.7) Version' of 'Linux 64-bit'.

### 1.2.3 RUN

Type './run\_ToTEM.sh <mcr\_directory>' at terminal to run it. For example, if you have version 9.9 of the MATLAB Runtime installed in the folder: '/mathworks/home/application/v99', run the shell by using this script,

```
./run_ToTEM_v3.sh /mathworks/home/application/v99
```

  
A blank space



# III. DETAILS

## ➤0. GUI

### ➤1. Import a crystal

### ➤2. Slicing dialog of multislice method

### ➤3. Parameter settings about simulation

3.1 Vibration, DW factors and absorption

3.2 Imaging region

3.3 Coherence

3.4 Lens parameters

3.5 Detector's angle format

## ➤4. Other functions

4.1 CPU/GPU parallels

4.2 HRTEM related

4.3 Save the file

## ➤5. Examples

### 5.1 HRTEM

5.1.1 HRTEM Simulation for one single image

5.1.2 Batch Simulation

### 5.2 CBED

### 5.3 STEM

### 5.4 IDPC

### 5.5 IDPC Simulation for MFI-type zeolites

# 0. GUI

ToTEM\_submit\_v1

Crystal HRTEM STEM/CBED Display Setting Contact

2. Crystal

Filename

Output  **2. Slicing**

Projected Potential

Potential param.

☐ Vibration

☐ Temperature Fa

☐ Lobato

☐ Height Influence in Slice

Imaging Region

Edge Size  Pixels **4. Show**

Imaging Resolution  A/Pixel

Imaging Region (Top-Left)

Angstrom

Width & Height of Imaging Region

Pixels

Scan Step  X Imaging Resolution

3. Mode

☐ HRTEM ☐ CBED ☐ CPU

☒ STEM ☐ IDPC ☒ GPU

Imaging Condition

Voltage:  kV Convergence:  mrad

Focus Spread:  nm Number:

Spatial Convergence:  mrad

Beam tilt:  mrad  degree

Magnitude(nm) & Angle(degree)

Defocus C1:  nm

2-fold astigmatism A1:  nm

3-fold astigmatism A2:  nm

Axial coma B2:  nm

(C3) Cs:  um

4-fold astigmatism A3:  um

Axial star aberration S3:  um

5-fold astigmatism A4:  um

3-lobe aberration D4:  um

coma B4:  um

C5:  mm

6-fold astigmatism A5:  mm

R5:  mm

S5:  mm

STEM Detector mrad

**5. Simulation**

Path  Save as:

Red: Imaging Region  
Black: Probe Size  
Green: Calculate Projected Potential

Angstrom

X-Extend


Y-Extend

Electron Incident (az & el degree)   **1. Top-View along Beams** **Top-view as screen**




# 1. Import a crystal

Import crystal information from a PDB file or an excel or text file.



Input Debye-waller (DW) factor for every element. The dialog box will pop up after importing a PDB file



**A \*.xls or \*.txt file will be generated after reading a pdb file.  
Next time you can import it without inputting DW factors.**


# Information in a xls or txt file.

	Atomic number	Coordinators (units: Å)			Valence state	Ionicity	temperature	Occupancy rate
	A	B	C	D	E	F	G	H
1	8	-1.973	0	-1.973	0	0	0.7875	1
2	8	-5.918	0	-1.973	0	0	0.7875	1
3	8	-9.863	0	-1.973	0	0	0.7875	1
4	8	-13.808	0	-1.973	0	0	0.7875	1
5	8	-17.753	0	-1.973	0	0	0.7875	1
6	8	-21.698	0	-1.973	0	0	0.7875	1
7	8	-25.643	0	-1.973	0	0	0.7875	1
8	8	-29.588	0	-1.973	0	0	0.7875	1

**You can conveniently edit it via other tools.**



After importing a model, click this button to view and calculate the model from the direction of electron incidence (as values in the edit boxes).




Or you can see the model from any angle after mouse rotating.

**Note: at least one button must be clicked**


## 2.Slicing dialog of multislice method

If some atoms were lost after slicing, please add a very small vacuum on top or bottom, such as 0.00001.



Click to show the slicing interface.

Necessarily!  
Click it to generate a new crystal



After setting the slicing parameters, you must click the button to display the slicing results.

Other format like:  
2,3,5(1),3(5:7)

To pick some slices from all slices


2	3	1	1	1	1	1	5	6	7
		5	6	7	5	6	7		

The content between brackets will be repeated several times.

**‘Top’ and ‘Bottom’ should be zero if repeating some slices!**



# Slicing result




# 3. Parameter settings about simulation

## 3.1 Vibration, Temperature Factors

$$B = \frac{8}{3} \pi^2 \langle dr^2 \rangle$$


$$V(\vec{r}) = \begin{cases} \int f(\vec{s}) \exp(-Bs^2) \exp(4\pi i \vec{s} \cdot \vec{r}) d^3s & \text{(iia)} \\ \int f(\vec{s}) \exp(4\pi i \vec{s} \cdot d\vec{r}) \exp(4\pi i \vec{s} \cdot \vec{r}) d^3s & \text{(iib)} \end{cases}$$

$f(\vec{s})$  is the parameterized scattering factor ignoring the TDS effect



(i)



- 
- **Vibration**: All atoms will vibrate according to the DW factor (Eq. i) and input vibration times in the frozen lattice mode. That is, use Eq. iib to simulate TDS.
  - **Temperature Factor and Frozen Lattice mode**: use Eq. iia.
  - **Height Influence in Slice**: Potential of each atom will contribute to multiple slices and  $\mu$  value.





☐ Height Influence in Slice 2

Projected Potential  
Potential param. Peng 0~6 A-1  
○ Vibration 30  
○ Temperature Fa  
Peng extended  
Lobato  
☐ Height Influence in Slice 2  
Imaging Region

Peng model and Peng extended model can take 'Height Influence in Slice', but the potential cannot be assigned in 3D domain if using Lobato model.

## 3.2 Imaging region

### STEM/IDPC Model

Imaging Region

Size  Pixels

Imaging Resolution  A/Pixel

Imaging Region (Top-Left)

Angstrom

Width & Height of Imaging Region


Pixels

Scan Step  X Imaging Resolution

3. Mode

☐ HRTEM ☐ CBED ☒ CPU

☒ STEM ☐ IDPC ☐ GPU



- Green box: region of the calculated potential, determined by the red and black regions
- Black box: probe wave's size is in 'Edge Size' and its physical size should be multiplied by 'Imaging Resolution'.
- 'Imaging resolution': the sampling rate of the probe wave.
- Red box: the imaging region by editing 'Imaging Region (Top-Left)' and 'Width & Height of imaging Region'.
- Scan Step: it multiplied by the sampling rate is the sampling rate in the red region.

# HRTEM Model

Imaging Region

Edge Size  Pixels

Imaging Resolution  Å/Pixel

Imaging Region (Top-Left)

Å


Width & Height of Imaging Region

Pixels

3. Mode

☒ HRTEM ☐ CBED ☐ CPU

☐ STEM ☐ IDPC ☐ GPU



- Green box: region of the calculated potential, determined by the red region and its edge, e.g. the side length is  $(512+64*2)$
- ‘Imaging resolution’: the sampling rate of the image.
- Red box: the imaging region by editing Imaging Region (Top-Left) and Width & Height of imaging Region. Its physical size should be multiplied by ‘Imaging Resolution’.



# CBED Model

Imaging Region

Size  Pixels

Imaging Resolution  A/Pixel


Imaging Region (Top-Left)

Angstrom

3. Mode

☐ HRTEM ☒ CBED ☐ CPU

☐ STEM ☐ IDPC ☐ GPU



- Green box: region of the calculated potential, and it is extended by a constant (4 Angstrom) to eliminate the aliasing.
- Red box: the imaging region by editing Imaging Region (Top-Left) and Width & Height of imaging Region. Its physical size should be multiplied by 'Imaging Resolution'.

## 3.3 Coherence


Imaging Condition

Voltage:	200	kV	Aperture	10	1/nm
Focus Spread:	3.0	nm	Number:	1	▼
Spatial Convergence:	0.3	mrad			

### HRTEM Model

- **Voltage**: Accelerated voltage
- **Temporal coherence**: represented by the root mean square of focal spread  $\Delta\varepsilon$  ('Focus Spread')
- **Spatial coherence**: represented by the semi-convergence angle  $\alpha$  of the incident electron beam ('Spatial Convergence')
- **Aperture**: Aperture of the objective lens

$$f_{\Delta}(m) = \frac{1}{\sqrt{2\pi}\Delta\varepsilon} \exp\left[-\frac{(m\delta\varepsilon)^2}{2\Delta\varepsilon^2}\right] \quad m \in [-M', \dots, 0, 1, \dots, M'] \text{ with } \sum_{m=-M'}^{M'} f_{\Delta}(m) = 1$$



The  $2*M+1$  is input in 'Number'.

Imaging Condition


Voltage: 200 kV Convergence 15 mrad

Focus Spread: 3.0 nm Number: 1

Spatial 0.3 mrad

## STEM/IDPC/CBED Model

- **The temporal coherence**: represented by the root mean square of focal spread  $\Delta\epsilon$  ('Focus Spread')
- **The semi-convergent angle**: the semi-convergence angle  $\alpha$  of the incident electron beam ('Convergence')
- **Voltage**: Accelerated voltage

$$f_{\Delta}(m) = \frac{1}{\sqrt{2\pi}\Delta\epsilon} \exp\left[-\frac{(m\delta\epsilon)^2}{2\Delta\epsilon^2}\right] \quad m \in [-M', \dots, 0, 1, \dots, M'] \text{ with } \sum_{m=-M'}^{M'} f_{\Delta}(m) = 1$$


The  $2*M+1$  is input in 'Number'.



## 3.4 Lens parameters

$$\xi = A_\xi \exp(-i\theta_\xi)$$


where  $\xi$  represents any quantity of astigmatism,  
and the wave aberration function  $\chi$  of lens is written as:

$$\begin{aligned} \chi(k, \Delta f, A_1, A_2, B_2, \dots) = & \operatorname{Re} \left\{ \frac{1}{2} \Delta f \lambda k k^* + \frac{1}{2} A_1 \lambda k^{*2} + \right. \\ & \frac{1}{3} A_2 k^2 k^{*3} + B_2 \lambda^2 k^2 k^* + \frac{1}{4} C_s \lambda^3 (k k^*)^2 + \frac{1}{4} A_3 \lambda^3 k^{*4} + \\ & S_3 \lambda^3 k^3 k^* + \frac{1}{5} A_4 \lambda^4 k^{*5} + D_4 \lambda^4 k^4 k^* + B_4 \lambda^4 k^3 k^{*2} + \\ & \left. \frac{1}{6} C_5 \lambda^5 (k k^*)^3 + \frac{1}{6} A_5 \lambda^5 k^{*6} + R_5 \lambda^5 k^5 k^* + S_5 \lambda^5 k^4 k^{*2} \right\} = \\ & \frac{1}{2} \Delta f \lambda k k^* + \frac{1}{4} C_s \lambda^3 (k k^*)^2 + \chi'(k, A_1, A_2, B_2, \dots) \end{aligned}$$

	$A_\xi$	$\theta_\xi$
Defocus C1:	0 nm	
2-fold astigmatism A1:	0 nm	0
3-fold astigmatism A2:	0 nm	0
Axial coma B2:	0 nm	0
(C3) Cs:	5 um	
4-fold astigmatism A3:	0 um	0
Axial star aberration S3:	0 um	0
5-fold astigmatism A4:	0 um	0
3-lobe aberration D4:	0 um	0
coma B4:	0 um	0
C5:	0 mm	
6-fold astigmatism A5:	0 mm	0
R5:	0 mm	0
S5:	0 mm	0


## 3.5 Detectors

➤ **STEM modes:** The Angle of the detector is separated by ‘,’ , 6 is inner angle and 20 is exterior angle. And different detectors is separated by ‘;’.



A screenshot of a software window titled "STEM Detector mrad". It features a yellow tab labeled "5. Simulation". Below the tab is a text input field containing the string "6,20;40,200;". To the right of the input field are up and down arrow buttons. At the bottom of the window, there is a "Path" button, a "Save as:" label, and a text field containing the word "test".


➤ **IDPC modes:** The Angle of the detector is separated by ‘,’ , 6 is inner angle and 20 is exterior angle. And different detectors is separated by ‘;’.



A screenshot of a software window titled "STEM Detector mrad". It features a yellow tab labeled "5. Simulation". Below the tab is a text input field containing the string "6,20;40,200;". To the right of the input field are up and down arrow buttons. At the bottom of the window, there is a "Path" button, a "Save as:" label, and a text field containing the word "test".

## 4. Other functions

### 4.1 STEM/IDPC/CBED related



➤ **HRTEM/CBED /STEM/IDPC modes**: users can use GPU/CPU multi-thread parallels.







2. Crystal

Filename

Output

In this input box enter the serial number of intermediate slices you want to output.



Note: The serial number must be an increasing sequence, and numbers are separated by ‘,’.




Phase of probe wave

Intensity profile of Probe Wave


Intensity of probe wave



Result




Blurred image




The STEM image can be blurred by a Gaussian function with a specified full width at half maxima (FWHM) and increasing the sampling, which is designed with reference to QSTEM.




## 4.2 HRTEM related





Phase plate





imag(CTF)



Real(CTF)




Damping(CTF)



## 4.3 Save the file

Select the  
folder to store  
the results



Name of  
result

The result will be named as **test\_STEM\_101.5037slice\_Dect\_40\_200.dat**


↓   ↓   ↓   ↓  
name   model   thickness   Detector angle

If the STEM image is blurred by a Gaussian function, the result will be named as **test\_STEM\_101.5037slice\_Dect\_40\_200\_source\_2\_zoom2.dat**

↓   ↓  
Value of source   Value of sampling

Open the files with the .dat suffix in matlab

```
fid=fopen('file name', 'r');  
data=fread(fid,[ size, size],'float');  
fclose(fid);  
data=data';
```




```
discontent -  
'Mid result _W20 * H20_Float Data'  
Mid result _W20 * H20_Float Data inC:  
|
```

Size: Output at the end of the program



## 5.Example



The simulation uses the  $\text{SrTiO}_3$ , the pictures is the slicing interface parameter setting, The slicing interface settings of all modes are the same. The interface parameter settings of different modes are different, and they will be given in different modes. The parameters of the potential calculation use Peng's parameters.

# 5.1 HRTEM

Crystal HRTEM STEM/CBED Display Setting

2. Crystal  
Filename: mp-5229\_SrTiO3101025.pc  
Output: 1 **2. Slicing**

Projected Potential  
Potential param: Peng 0~6 A-1  
☐ Vibration 1  
☐ Temperature Factor  
☐ Height Influence in Slice 2

Imaging Region  
Size: 64 Pixels **4. Show**  
Imaging Resolution: 0.1 A/Pixel  
Imaging Region (Top-Left): -3 -3 Angstrom

3. Mode  
☐ HRTEM ☒ CBED  
☐ STEM ☐ IDPC  
☒ CPU ☐ GPU


Imaging Condition  
Voltage: 200 kV Convergence: 10 mrad  
Focus Spread: 3.0 nm Number: 1  
Spatial Convergence: 0.3 mrad  
Beam tilt: 0 mrad 0 degree  
Magnitude(nm) & Angle(degree)  
Defocus C1: 0 nm  
2-fold astigmatism A1: 0 nm 0  
3-fold astigmatism A2: 0 nm 0  
Axial coma B2: 0 nm 0  
(C3) Cs: 20 **um**  
4-fold astigmatism A3: 0 **um** 0  
Axial star aberration S3: 0 **um** 0  
5-fold astigmatism A4: 0 **um** 0  
3-lobe aberration D4: 0 **um** 0  
coma B4: 0 **um** 0  
C5: 0 **mm**  
6-fold astigmatism A5: 0 **mm** 0  
R5: 0 **mm** 0  
S5: 0 **mm** 0

**5. Simulation**  
Path Save as: test

→ interface



Thickness is 101.5037Å

result




## 5.1.2 Batch HRTEM Simulation

**Application:** The simulation image can be produced in large quantities within a certain astigmatic range



Exit wave function  
generated in HRTEM mode






Please do not change the format of this file!



```
parameter_batchsimulation.txt - 记事本
文件(F) 编辑(E) 格式(O) 查看(V) 帮助(H)
4 % myrandrange.focus=allvalue(1) %units: nm will be chanced with both positive and nagative directions
2.5 %myrandrange.A1=allvalue(2) %units:nm
180 %myrandrange.phiA1=allvalue(3); %units:degrees max 180
10 %myrandrange.A2=allvalue(4); %units:nm
180 %myrandrange.phiA2=allvalue(5); %units:degrees max 180
10 %myrandrange.B2=allvalue(6); %units:nm
180 %myrandrange.phiB2=allvalue(7); %units:degrees max 180
1000 %myrandrange.Cs=allvalue(8); %unit nm ****
700 %myrandrange.A3=allvalue(9); %units:nm ****
180 %myrandrange.phiA3=allvalue(10); %units:degrees max 180
600 %myrandrange.S3=allvalue(11); %units:nm ****
180 %myrandrange.phiS3=allvalue(12); %units:degrees max 180
50000 %myrandrange.A4=allvalue(13); %units:nm ****
180 %myrandrange.phiA4=allvalue(14); %units:degrees max 180
0 %myrandrange.D4=allvalue(15); %units:nm ****
0 %myrandrange.phiD4=allvalue(16); %units:degrees max 180
0 %myrandrange.B4=allvalue(17); %units:nm ****
0 %myrandrange.phiB4=allvalue(18); %units:degrees max 180
0 %myrandrange.A5=allvalue(19); %units:nm ****
0 %myrandrange.phiA5=allvalue(20); %units:degrees max 180
0 %myrandrange.C5=allvalue(21); %units:nm **** %其实也可以设置百分比。
```



Import ranges to randomly generate values around each aberration.



Beam tilt:	0	mmrad	0	degree
Magnitude(nm) & Angle(degree)				
Defocus C1:	-2.5	nm		
2-fold astigmatism A1:	2.5	nm	0	
3-fold astigmatism A2:	40	nm	0	
Axial coma B2:	35	nm	0	
(C3) Cs:	4	um		
4-fold astigmatism A3:	0.8	um	0	
Axial star aberration S3:	0.8	um	0	
5-fold astigmatism A4:	50	um	-130.9	
3-lobe aberration D4:	32.4	um	112.1	
coma B4:	60	um	0	
C5:	2.71	mm		
6-fold astigmatism A5:	2.562	mm	-65.9	
R5:	0	mm	0	
S5:	0	mm	0	

**Image's number of:** The number of images generated

**Process number:** Number of processes

**Batch size for one process:** The number of images generated per batch per process

Input Ima... — □ ×

Enter wave size: width (Cubic)  
320

Input flag: Show wave(S) and Run simulation(R)  
R

Images' number of  
500

Process number  
2

Batch size for one process  
5


TOP-LEFT coordinate (x-y)  
101 51

Intercept size (width & hight)  
50 70


确定 取消

**Run simulation(R):** start to generate images

**Show wave(s):** display the phase of wave and show the interception area



# Batch simulation result



**Fifty images of HRTEM model were successfully generated!**



## 5.2 CBED

Crystal HRTEM STEM/CBED Display Setting

2. Crystal  
Filename: mp-5229\_SrTiO3101025.pc  
Output: 1 **2. Slicing**

Projected Potential  
Potential param.: Peng 0~6 A-1  
☐ Vibration 1  
☐ Temperature Factor  
☐ Height Influence in Slice 2

Imaging Region  
Size: 256 Pixels **4. Show**  
Imaging Resolution: 0.1 A/Pixel  
Imaging Region (Top-Left): 10 10 Angstrom

3. Mode  
☐ HRTEM ☒ CBED  
☐ STEM ☐ IDPC  
☒ CPU ☐ GPU


Imaging Condition  
Voltage: 300 kV Convergence: 12 mrad  
Focus Spread: 3.0 nm Number: 1  
Spatial Convergence: 0.3 mrad  
Beam tilt: 0 mrad 0 degree  
Magnitude(nm) & Angle(degree)  
Defocus C1: 0 nm  
2-fold astigmatism A1: 0 nm 0  
3-fold astigmatism A2: 0 nm 0  
Axial coma B2: 0 nm 0  
(C3) Cs: 20 um  
4-fold astigmatism A3: 0 um 0  
Axial star aberration S3: 0 um 0  
5-fold astigmatism A4: 0 um 0  
3-lobe aberration D4: 0 um 0  
coma B4: 0 um 0  
C5: 0 mm  
6-fold astigmatism A5: 0 mm 0  
R5: 0 mm 0  
S5: 0 mm 0

**5. Simulation**  
Path Save as: test




interface

result






## 5.3 STEM




Thickness is 101.5037Å


→ interface




result



blurred by Gaussian Function




## 5.4 IDPC




Thickness is 23.7Å

→ interface


result



blurred by Gaussian Function



## 5.5 IDPC Simulation for MFI-type zeolites




Thickness is 38.3Å

← interface


Structure

010 direction






result



blurred by Gaussian Function



Prof. Fang Lin developed the codes on CPU and GUI. Mr. PengJun Yuan at Xiangtan University developed the CUDA codes to speed up calculation, and Mr. KePeng Wu tested the program and developed some codes on CPU. Prof. Yuan Yao and Prof. ChuanHong Jin offers theoretical support. Mr. ZhiQun Li at South China University of Technology have done some previous work.

We thank our collaborators for their suggesting: Prof. Yu Han at KAUST, Prof. JunHao Lin at Southern University of Science and Technology, and Prof. DaLiang Zhang at Chongqing University.

We acknowledge the financial support from the National Science Foundation of China.



**Thank you for using ToTEM**